

Epinephrine Infusions

| Saline Bag | 250 mL | 500 mL | 1000 mL |
|------------|--------------------------|--------------------------|--------------------------|
| Drip Set | 60 gtts | 10 gtts | 10 gtts |
| Add Drug | 1.0 mg Epinephrine | 0.5 mg Epinephrine | 1.0 mg Epinephrine |
| Result | 4 mcg/mL | 1 mcg/mL | 1 mcg/mL |
| 1 mcg/min | 15 gtts/min (1 gt/4s) | 10 gtts/min (1 gt/6s) | 10 gtts/min (1 gt/6s) |
| 2 mcg/min | 30 gtts/min (1 gt/2s) | 20 gtts/min (1 gt/3s) | 20 gtts/min (1 gt/3s) |
| 4 mcg/min | 60 gtts/min (1 gt/1s) | | |
| 6 mcg/min | | 60 gtts/min (1 gt/1s) | 60 gtts/min (1 gt/1s) |
| 8 mcg/min | 120 gt/min (2 gt/1s) | | |
| 12 mcg/min | | 120 gt/min (2 gt/1s) | 120 gt/min (2 gt/1s) |

